

Technical Guides

Environmental Test Chambers & Associated Systems

ENGINEERING REFERENCE

Comprehensive engineering specifications, calibration metrics, and mechanical system layout instructions for environmental test chambers, refrigeration systems, humidity systems, airflow circuits, controls, utilities and connected service infrastructure.

Designed for service engineers, maintenance teams, facility engineers, calibration personnel and equipment owners.

Scope

- Temperature and humidity chamber engineering fundamentals
- Refrigeration and air-circulation system architecture
- Electrical, control and sensor reference points
- Calibration methodology and performance metrics
- Mechanical installation and utility-layout guidance
- Commissioning, diagnostics and documentation requirements

IMPORTANT

Values in this guide are practical engineering reference ranges, not substitutes for the chamber manufacturer's approved drawings, component data, safety instructions or the applicable test standard. Final settings and acceptance criteria must be verified against the specific equipment and application.

1. Chamber System Architecture

An environmental test chamber is a controlled thermal system in which refrigeration, heating, humidity, airflow, sensing and control must operate as one integrated loop. Reliable performance depends as much on mechanical layout and airflow integrity as on individual component capacity.

1.1 Functional system blocks

System	Primary Function	Typical Elements
Test space	Controlled environment around DUT	Insulated enclosure, door, viewing window, ports, racks
Air circulation	Uniformly distributes conditioned air	Blower, motor, wheel, baffles, supply/return path
Heating	Raises chamber temperature	Electric heaters, SSR/contactor, over-temperature protection
Refrigeration	Removes sensible/latent heat	Compressor, condenser, expansion device, evaporator, controls
Humidity	Adds/removes moisture	Steam/water source, dehumidification coil, drain, humidity sensor
Controls	Regulates setpoint and sequences equipment	PLC/controller, HMI, I/O, relays, safety interlocks
Measurement	Provides control and verification data	RTD/thermocouple, RH sensor, pressure transducers
Utilities	Supports continuous operation	Power, cooling water/air, DI water, drain, network

1.2 Airflow principle

- Maintain a clear supply-to-return path. Avoid placing the test specimen directly against supply or return openings.
- Do not overload shelves or block perforations; restricted airflow increases gradients and recovery time.
- Sensor placement should represent the controlled air condition, not direct heater, evaporator or specimen radiation.
- After any blower, baffle or duct modification, re-check spatial uniformity before returning the chamber to validated service.

FIELD RULE

A chamber that reaches setpoint can still be out of specification. Always distinguish controller indication from actual spatial uniformity, fluctuation and recovery performance.

2. Engineering Specification Reference

Use this section to define the technical baseline for an installed chamber. Record the actual OEM values where available; use “N/A” rather than assumptions.

Parameter	Engineering Reference	Record / Verify
Temperature range	Application-specific; verify rated minimum and maximum	____ °C to ____ °C
Humidity range	Defined only within manufacturer temperature/RH envelope	____ %RH to ____ %RH
Temperature fluctuation	Short-term variation at a stable point	____ °C
Temperature uniformity	Spatial difference across defined measurement points	____ °C
Ramp rate	Heating/cooling rate under defined load and range	____ °C/min
Recovery time	Time to return within tolerance after disturbance	____ min
Workspace volume	Usable internal test volume	____ L / m ³
Electrical supply	Voltage, phase, frequency, connected load	_____
Cooling method	Air-cooled / water-cooled / remote condenser	_____
Refrigerant circuit	Single-stage / cascade / other	_____
Controller	Make, model, firmware/version	_____
Communications	Ethernet / RS-485 / protocol / remote monitoring	_____

2.1 Performance terminology

Metric	What it tells you
Fluctuation / stability	How much temperature or humidity varies with time at a fixed point after stabilization.
Uniformity / gradient	How evenly the environment is distributed across the usable workspace.
Deviation / error	Difference between chamber indication and a traceable reference measurement.
Ramp rate	Rate of temperature change; specify whether average, linear-zone, loaded or unloaded.
Recovery	Ability to return to tolerance after door opening, load insertion or programmed transition.

3. Calibration & Measurement Guide

Calibration should establish traceability and quantify measurement performance under defined conditions. The procedure, sensor count, locations, dwell time and acceptance limits must match the applicable standard, customer requirement and chamber use case.

3.1 Temperature mapping

- Use calibrated reference sensors with uncertainty appropriate to the required chamber tolerance.
- For spatial mapping, distribute sensors through the defined working volume; a common engineering approach uses corners plus a centre point, with additional points for large chambers.
- Keep sensors away from walls, direct supply jets, heaters and evaporator surfaces unless the procedure specifically requires those locations.
- Allow stabilization before collecting data. Record sufficient samples to characterize both mean value and short-term variation.

Metric	Calculation / Interpretation	Report
Indication error	Chamber display minus reference value at control location	Error at each setpoint
Spatial uniformity	Difference between highest and lowest mapped location values over defined period	Worst-case spread
Temporal fluctuation	Variation at an individual sensor after stabilization	Peak-to-peak or defined method
Measurement uncertainty	Combined uncertainty of reference system and method	Expanded uncertainty + coverage factor
Correction	Applied or documented offset where permitted	Before/after values

3.2 Humidity calibration

- Verify humidity only inside the chamber's specified controllable temperature/RH envelope.
- Use a traceable reference hygrometer or dew-point instrument suitable for the range.
- Allow longer stabilization at low humidity and after large humidity transitions.
- Inspect water quality, steam generation, drains and wet surfaces before attributing RH error solely to the sensor.

CALIBRATION RECORD

Record setpoint, chamber indication, reference value, correction/error, stability period, sensor serial number, calibration due date, uncertainty, ambient conditions and pass/fail criteria.

4. Refrigeration System Guide

The refrigeration circuit determines low-temperature capability, pull-down performance and heat-removal capacity. Service decisions should be based on measured pressures, temperatures and component condition rather than pressure alone.

4.1 Typical single-stage flow

COMPRESSOR → CONDENSER → RECEIVER / FILTER-DRIER → EXPANSION DEVICE → EVAPORATOR → SUCTION

4.2 Cascade system concept

Low-temperature chambers may use two refrigeration stages coupled through a cascade heat exchanger. The high stage rejects heat from the low stage; stable operation therefore depends on correct sequencing and adequate high-stage heat rejection before the low stage is heavily loaded.

Check	Engineering Observation
Compressor	Current, oil condition/level where visible, discharge temperature, noise/vibration
Condenser	Cleanliness, fan/pump operation, entering/leaving air or water condition
Expansion device	Stable feed, superheat response, icing pattern, sensing bulb/contact where applicable
Evaporator	Even frost/temperature pattern, airflow, defrost effectiveness, drain condition
Filter-drier / sight glass	Restriction indicators, moisture indication where fitted, abnormal temperature drop
Piping	Oil return, insulation, vibration support, abrasion, leaks, service-valve condition

4.3 Diagnostic measurements

- Suction and discharge pressure
- Suction-line temperature and calculated superheat
- Liquid-line temperature and calculated subcooling where applicable
- Compressor current versus nameplate/expected operating condition
- Condenser entering/leaving temperature
- Evaporator air-on/air-off temperature
- Cascade interstage temperature/pressure for cascade systems

SAFETY

Refrigerant handling, electrical testing and pressure-system work must be performed by qualified personnel using the correct PPE, recovery equipment and procedures for the refrigerant and system design.

5. Humidity, Airflow & Water Systems

5.1 Humidity generation and removal

Humidity control is a balance between moisture addition, cooling/dehumidification, chamber leakage and test-load moisture effects. Poor water quality or drainage can create scale, contamination, unstable RH and sensor drift.

Area	Technical Requirement
Water supply	Use quality specified by OEM; DI/demineralized water is commonly required for steam/humidity systems.
Steam generator	Inspect heater surface, level control, scale and electrical connections.
Drain	Maintain free drainage, trap arrangement and temperature-compatible piping.
Dehumidification	Verify refrigeration/dehumidification coil operation and condensate removal.
RH sensor	Keep clean, protected from direct spray/condensation and within rated temperature range.

5.2 Air-circulation mechanical layout

- Blower wheel rotation must match design direction.
- Motor current and bearing condition should be checked after mechanical work.
- Baffles and plenums must be refitted in their original orientation and sealed where required.
- Maintain clearance around return air and supply air paths.
- Large test specimens may require a defined minimum clearance on all sides to prevent recirculation pockets.

5.3 Typical airflow verification

Observation	Possible Cause
Large top-to-bottom gradient	Blocked return, weak circulation, specimen obstruction, baffle issue
Slow recovery	Insufficient airflow, refrigeration/heating capacity issue, excessive load
Local hot/cold zone	Direct jet, poor baffle distribution, sensor placement, specimen heat load
RH instability	Air leakage, water/steam issue, refrigeration cycling, sensor contamination

6. Electrical, Controls & Safety Systems

6.1 Electrical inspection reference

Item	Verify
Incoming power	Voltage, phase balance, frequency, protective earth, breaker rating
Panel	Cleanliness, cooling fans/filters, terminal condition, signs of overheating
Heater circuit	Contactor/SSR operation, current balance, insulation condition, over-temp cutout
Motors	Current, overload setting, rotation, VFD status where fitted
Safety chain	Emergency stop, door interlock, pressure switches, thermal protection, alarms
Control power	Power supply output, PLC/controller status, I/O indication, relay operation
Communications	Ethernet/serial link, remote alarms, time synchronization, data logging

6.2 Sensor and control loop checks

- Compare controller process value with an independent calibrated reference.
- Inspect sensor wiring, terminals, shielding and routing near high-current conductors.
- Confirm configured sensor type and scaling after controller, PLC or transmitter replacement.
- Record controller offsets before changing them; calibration correction should be traceable.
- Verify independent safety cut-outs separately from the normal control loop.

DO NOT BYPASS SAFETIES

Pressure switches, thermal cut-outs, overloads, emergency stops and interlocks are protective devices. Temporary bypasses should not be used to return equipment to normal service.

7. Mechanical Layout & Installation Instructions

A good installation provides service access, stable utilities, drainage, ventilation and a clear path for future component replacement. Layout should be reviewed before chamber placement, not after commissioning.

7.1 Site-layout checklist

Layout Item	Design / Verification Requirement
Equipment clearance	Provide OEM-required service and ventilation clearance around chamber and condensing units.
Floor / foundation	Level, load-bearing surface; verify anchoring or anti-vibration provisions where required.
Door access	Allow full door swing and safe movement of test specimens.
Remote refrigeration	Route piping to control pressure drop, oil return, vibration and thermal expansion.
Cooling water	Provide isolation, strainer, flow/pressure suitability and accessible connections.
Drain	Continuous fall where gravity drained; avoid back-pressure and unintended traps.
Humidity water	Dedicated clean supply with accessible isolation and filtration/treatment as specified.
Electrical	Accessible isolator, correct cable sizing/protection, earthing and phase sequence.
Exhaust / heat rejection	Prevent recirculation of condenser discharge air and excessive plant-room temperature.
Network / IoT	Provide reliable Ethernet/Wi-Fi/cellular path as applicable and agreed IT access.

7.2 Serviceability rules

- Keep compressors, valves, filters, electrical panels and sensors accessible without dismantling unrelated building work.
- Use pipe supports that prevent vibration transfer and rubbing.
- Insulate cold suction lines and other surfaces where condensation is expected.
- Label isolation points, flow direction, circuits and remote equipment.
- Leave space for future compressor, motor, heater and controller replacement.

8. Commissioning & Acceptance

8.1 Pre-start verification

- Confirm chamber is level, secure and free of transport restraints.
- Verify electrical supply, earthing, phase sequence and protective devices.
- Verify condenser airflow or cooling-water flow.
- Check refrigeration piping and joints for signs of leakage or damage.
- Confirm humidity water and drain arrangements.
- Inspect blower rotation, baffles, ports, door seals and safety devices.
- Confirm controller configuration, sensor inputs and programmed limits.

8.2 Functional acceptance sequence

Stage	Test	Record
1	Power-up and safety-chain verification	Alarms / interlocks / E-stop
2	Ambient-temperature circulation check	Blower current, noise, airflow
3	Heating test	Ramp, heater current, over-temp function
4	Cooling test	Pull-down, pressures/temperatures, compressor current
5	Humidity test where applicable	RH generation, dehumidification, drain
6	Stability test	Temperature/RH fluctuation after stabilization
7	Spatial mapping if required	Uniformity across defined workspace
8	Communications / logging	Data capture, alarms, remote connection
9	Documentation handover	Drawings, manuals, calibration, service baseline

8.3 Baseline data to preserve

- No-load pull-down and heat-up curves
- Normal compressor current and key operating pressures/temperatures
- Temperature and RH calibration results
- Controller configuration and firmware version
- PLC/HMI backup where permitted
- Electrical and refrigeration drawings
- Installed component make/model/serial numbers
- Photographs of panel, piping and sensor locations

9. Diagnostic Guide

Symptom	Check First	Then Investigate
Does not cool	Setpoint, alarms, compressor enable, condenser condition	Refrigeration pressures/temps, expansion device, charge/restriction
Slow pull-down	Load, airflow, condenser, door seal	Capacity, evaporator feed, compressor condition, insulation
Temperature unstable	Airflow, sensor, control output	PID/control settings, SSR/contactors, refrigeration cycling
High temperature gradient	Specimen blockage, blower, baffles	Air distribution, sensor locations, heat load
Humidity low	Water supply, generator, leaks	Heater/level control, sensor, refrigeration/dehumidification
Humidity high	Drain, dehumidification, sensor wetting	Cooling coil, control logic, leakage/condensation
High discharge pressure	Condenser cleanliness/airflow/water flow	Overcharge, non-condensables, ambient/load
Compressor trips	Supply, overload, condenser	Pressures, current, winding condition, oil/refrigerant return
Repeated sensor error	Connections, configured sensor type	Cable shielding, transmitter, controller input

DIAGNOSTIC PRINCIPLE

Capture the condition before changing it: setpoint, process value, alarm history, currents, pressures, temperatures and recent service events. A good service record turns repeat faults into diagnosable patterns.

10. Technical Service Record

Use this page as a minimum technical record for preventive maintenance, troubleshooting, calibration or retrofit work.

Equipment Detail	Record
Customer / site	
Equipment / asset ID	
Manufacturer / model	
Serial number	
Temperature range	
Humidity range	
Controller / firmware	
Refrigeration type / refrigerant	
Date / engineer	

Measured operating data

Parameter	Before Service	After Service
Chamber temperature	____ °C	____ °C
Chamber RH	____ %RH	____ %RH
Suction pressure	____	____
Discharge pressure	____	____
Suction line temperature	____ °C	____ °C
Liquid line temperature	____ °C	____ °C
Compressor current	____ A	____ A
Blower current	____ A	____ A
Supply voltage	____ V	____ V

Findings & recommendations

Condition / fault observed

Corrective work completed

Parts replaced

Further action / quotation required

11. Levito Technical Support Framework

Technical guides are most useful when they form part of a complete equipment lifecycle record. Levito can use the same engineering baseline across service, preventive maintenance, calibration, modernisation and connected monitoring so that every intervention builds on the equipment's previous history.

Lifecycle Layer	Technical Output
Levito Service	Fault diagnosis, repair data, parts history, commissioning records
Levito Care	PM schedules, calibration status, compliance records, SLA history
Levito Modernise	Obsolescence review, controls/electrical/refrigeration upgrades, refurbishment
Levito Connect	Live parameters, alarms, uptime, trends, remote diagnostics and maintenance alerts

ONE ASSET. ONE TECHNICAL HISTORY.

Maintain a digital record for each chamber: specifications, drawings, calibration certificates, service history, replaced parts, operating trends, alarms and recommendations. This creates the foundation for condition-based and predictive maintenance.

For chamber-specific engineering support, calibration planning, retrofit assessment or connected monitoring, contact Levito through levito.in.